

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®


This document will review the updates and enhancements that have been made to Rental Products and Rental Processes in MotionMD to make the Rental Process more user friendly. This document contains the following sections:

1. Making an Existing Product Rentable
2. Adding a New Rental Product
3. Rental Products List Page Overview
4. Rental Product Page Overview
5. Dispensing a Rental Product Manually
6. Dispensing a Rental Product from a Patient Agreement
7. Dispensing a Rental Product from a Patient Agreement using a Protocol
8. Dispensing a Rental Product on a Patient Agreement from the Rental Products tab
9. Completing the Rental Cycle
10. Servicing a Rental Product
11. Inactivating a Rental Product
12. Additional Information


Prior to adding a New Rental Product, a Product and Inventory Item must be created. Follow the training document for creating new Products and Inventory Items. Once those steps are completed add Stock Records for each Inventory Location Rental Products will be dispensed from and set Par Level to Zero and On-hand to 1,000.

1. To Make an Existing Product Rentable

On the **Product** page, filter to **Account** and select the product.




Once on the Product select the **Rentals** tab.




OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

On **Rental Details**, click the **Edit** button




Check the **Rentable?** checkbox and click **Save Changes**. Checking Rentable allows this Product to be created as a Rental Product.




2. Adding a New Rental Product

From the **Inventory** tab, select the **Rental Products** option




OVERVIEW OF RENTAL PROCESS IN MOTIONMD®


This brings you to the **Rental Products list** page, to add a new Rental Product, click the **New Rental Product** button



Search for and Select the Product to be Rented.



Product Tooltip:



Enter the **Serial Number** of the Product and select the **Inventory Location** coinciding with a location the Stock Record was created for. These two fields are mandatory. Click the **Create Rental Product** button.

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

Rental Product New

Product 📘 Spectra Knee CPM ✕

Rental Unit Name

* Serial Number 📘 12345

Asset tag

Inventory location CSA Ortho (South) ▼

Notes

[Create Rental Product](#) [Cancel](#)

Serial Number Tooltip:

Rental Product New

Product 📘 Select a Product

Rental Unit Name

* Serial Number 📘

Asset tag

Serial Numbers on rental products cannot be updated. If this product is no longer in service, uncheck the Active? check box below and Save.

A message will show at the top of the page that the Rental Product was successfully created. Click the **Rental Products** link.

MOTIONMD Clinical Inventory Revenue Reports Admin Help

Spectra Knee CPM was successfully created.

[Rental Products](#) Spectra Knee CPM (12345)

Spectra Knee CPM Available [Create Patient Agreement](#) [View All Activity](#) [Edit Product](#)

This product is **Active** and currently **Available**

Product is not currently rented
Complete Service to add product to a PA or manually check out from this page [Service](#) [Check Out](#)

Rental History (0 Rentals / 0 Services)

Details


Asset Tag

Serial Number 12345

Product Spectra Knee CPM

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

The new newly created rental product will show up on the Rental Products List page and is Available for rental.



3. Rental Products List Page Overview


The **Rental Products List** page shows all Active Rental Products and has Search and Filter capabilities along with color coding notifications

Search box will work using full or partial Product Name or Serial Number


Status dropdown allows filtering the Rental list by any Status option

Active filter is defaulted and can be changed to filter to Inactive units

Inventory Locations filter is defaulted to All Inventory Locations. The dropdown allows filtering to a specific location



Color coding will highlight Rental Products that are Available, Currently Rented, Due Soon and Past Due. The Status of the Rental Product can be updated from this page by clicking the button on the row.



OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

4. Rental Product Page Overview

The **Rental Product** page shows the current and historical rental cycles for that particular unit (serial number).

The top section shows the current rental cycle and its status. Directly below is a list of all **Previous Rental Cycles** that will clearly indicate whether the unit was manually dispensed or dispensed from a Patient Agreement (see next two sections). The Status can be updated from this page by clicking the button in the current cycle section at the top.

The screenshot shows the 'Rental Products / Knee CPM (23456)' page. At the top, there are buttons for 'Create Patient Agreement', 'View All Activity', and 'Edit Product'. The main content area is divided into two sections: 'Currently Rented' and 'Rental History'. The 'Currently Rented' section shows the patient 'Bob Smith' (45 Day) and a 'Rental Cycle' that is 'Past Due'. It includes a table with columns for 'Checked Out', 'Picked Up', and 'Cleaned'. The 'Checked Out' row shows it was checked out by Leigha Wismann on 04/07/2022 at 04:48 PM. The 'Picked Up' row shows a warning icon and the text 'Expected to be returned on 05/22/2022 04:48 PM'. There is a 'Pick Up' button. The 'Rental History' section shows 'Manually Checked Out' (3 Day) on 11/29/21 by Kristen DiRocco. A 'Details' sidebar on the right lists 'Asset Tag', 'Serial Number', 'Product', 'Manufacturer', and 'Inventory Location'.

5. Manually Dispensing a Rental Product


Use this process if Patient Agreements in MotionMD are not going to be used as AOBs for the Rental Products.

From the **Rental Products list** page, find the product matching the serial number being dispensed and click the **Check Out** button or click into the **Rental Product** and click the **Check Out** button. A popup box will appear, enter the **number of days** the unit is being rented, this is a mandatory field, and click **Submit** button.


The screenshot shows the 'Rental Products' list page. A 'Check Out' popup box is open over the 'Knee CPM' product. The popup has a 'Rental period' field with a red box around it, containing '# of days'. Below the field is an 'Add Notes (Optional)' text area. At the bottom of the popup are 'Cancel' and 'Submit' buttons, with a red box around the 'Submit' button. In the background, the 'Rental Products' list is visible, showing 'Spectra Knee CPM' and 'Knee CPM' products. A 'Check Out' button is highlighted with a red box on the 'Knee CPM' product row.

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

The **Rental Product** will show as **Currently Rented** on the list page.



Within the **Rental Product** page, it shows the unit was manually checked out and the number of days of the rental period. The **Checked-Out** user and date and **Expected Pick Up** date are listed under the **Rental Cycle** section.




6. Dispensing a Rental Product from a Patient Agreement

To create a **Patient Agreement** for patient being dispensed a Rental Product ensure the Patient Agreement is created within the Clinic Location the Rental Product was created under.

On the **Products** section on the Patient Agreement the Rental Product can be searched by the **Product Name** or **Serial Number**.

If searching by **Name**, all Serial Numbers created under that product will be listed under a new **Rental Products** tab and the availability of each Serial Number is indicated.

If searching by **Serial Number**, only a product created matching that serial number will be displayed. Partial Serial Number searches will also bring back results.



OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

Dispensing a Rental Product from a Patient Agreement Cont.

Click the **Add Item(s)** button for the serial item being dispensed and complete the product form by adding **Involved Side**, **Rental Period Dates** (both mandatory fields) then click **Add Selected Item(s)**.

Note, the Serial Number field cannot be altered, if this is not the correct serial number, click Cancel and search for the correct serial number.

Add Spectra Knee CPM

Account
Spectra Knee CPM
DJO - Donjoy • Knee • CPM

Spectra Knee CPM (23456) Add Item(s)

Inventory Items 1 Selected Add Selected Item(s)

Select Options Required fields are marked with an asterisk *

98765 ✕
Available

Involved Side * Self Pay/TOS *

Left Right N/A No Yes

Fulfill From: Drop Ship

* Inventory Location: Ivory Ortho

* Billing Channel: CSA Ortho

Serial #: 98765

* Rental Period: 21 Start Rental on Date of Surgery

Cancel Add Selected Item(s)

The Product displayed on the Patient Agreement will indicate the Rental Period and the Serial Number dispensed. Clicking the View Rental Product link will take you to that Rental Products page.

Products Scan Barcode

(12345 / DJO - Donjoy)

Rental - 21 day 98765

[View Rental Product](#)


Self Pay Price: \$0.00 Billed Charge: \$0.00	Billing Channel: CSA Ortho Self Pay/TOS: No
---	--

Category
CPM
Side
Left
HCPCS


Edit Remove

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®


As soon as the **Patient signs** the Patient Agreement, the Rental Product will change to **Checked Out**.



Hovering over the tooltip will provide details on the rental product's current status, not necessarily its stats on the PA.



Navigating to that Rental Products page will show the unit as **Currently Rented/Checked Out** and the **Rental Cycle** will indicate the Patient it's currently rented to, the Rental Period in Days, PA number, MRN (if added) and Patient DOB.



OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

7. Dispensing a Rental Product on a Patient Agreement using Provider Protocol

Create a **Patient Agreement** for patient being dispensed a Rental Product using **Provider Protocol**. Ensure the Patient Agreement is created within the Clinic Location the Rental Product was created under.


Select the Provider from the drop-down menu and **Use Protocol > Choose Protocol > Select Serial Number**.

Once the PA is created, the Rental Product selected on the Protocol will auto-populate. The user will need to add **Involved Side**, **Rental Period Dates** (both mandatory fields) then click **Save Changes**.

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

8. Dispensing a Rental Product on a Patient Agreement from the Rental Products tab

On **Rental Products**, select the Rental Product and **Serial Number** for the Rental item dispensed. In the Rental Product select **Create Patient Agreement**.




Once the PA is created, the Rental Product selected on the Protocol will auto-populate. The user will need to add **Involved Side**, **Rental Period Dates** (both mandatory fields) then click **Save Changes**.

9. Completing the Rental Cycle Process

Whether the unit was Manually Checked Out or Checked Out from a Patient Agreement the remaining steps to complete the Rental Cycle are the same and can be completed from the List page or within the Rental Product page.

After a unit is **Checked Out**, the next status is **Pick Up**. When the unit has been picked up or returned, click the **ePicked Up** button.



OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

The next status is **Clean**, which indicates the unit is ready to be cleaned after use. You will also see the Picked-Up section has updated with the Date and User who selected the Picked Up button. When unit is ready to be cleaned, click the **Clean** button.

Currently Rented

Create Patient Agreement

Spectra Knee CPM (23456)

This product is **Active** and currently **Checked Out**

Manually Checked Out 1 Day
04/08/25 by Nicole Valentin

Add Note

Clean

Rental Cycle	
Checked Out	by Nicole Valentin on 04/08/2025 06:42 PM
Picked Up	by Nicole Valentin on 04/09/2025 09:14 AM
Cleaned	

The next status is **Check In**. You will now see the Date and User who clicked the Clean button. Once the unit has been cleaned and is ready to be rented, click the **Check In** button. Until the Check In button is clicked, this unit cannot be checked out for a new rental cycle or added to a Patient Agreement.

Currently Rented - Cleaning

Create Patient Agreement

Spectra Knee CPM (23456)

This product is **Active** and currently **Checked Out**

Manually Checked Out 1 Day
04/08/25 by Nicole Valentin

Add Note

Check In

Rental Cycle	
Checked Out	by Nicole Valentin on 04/08/2025 06:42 PM
Picked Up	by Nicole Valentin on 04/09/2025 09:14 AM
Cleaned	by Nicole Valentin on 04/09/2025 09:15 AM

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

After the **Check In** button is clicked, the unit will show as **Available** and ready for **Check Out**. The rental that was just Checked In will show at the top of the list under Previous Rental Cycles.

Spectra Knee CPM Available Create Patient Agreement

This product is **Active** and currently **Available**

Product is not currently rented
Complete Service to add product to a PA or manually check out from this page

Service Check Out

Rental History (2 Rentals / 0 Services)

Joe Test 14 Day		ID 107742
08/14/20 by Nicole Valentin		MRN -
		DOB 1957-06-11
Checked Out	by Nicole Valentin on 08/14/2020 06:14 PM	
Picked Up	by Nicole Valentin on 08/14/2020 06:19 PM	
Cleaned	by Nicole Valentin on 08/14/2020 06:21 PM	
Checked In	by Nicole Valentin on 08/14/2020 06:22 PM	

10. Servicing a Rental Product

If a Rental Product needs to be Serviced (repaired), the unit can be set to an in **Service** status on the Rental Product page. Note, the Service button will only be displayed if the unit is in the Available status. From the Rental Product page, click the **Service** button.

Spectra Knee CPM Available Create Patient Agreement

This product is **Active** and currently **Available**

Product is not currently rented
Complete Service to add product to a PA or manually check out from this page

Service Check Out

The status will update to Active and currently In Service. When servicing has been completed, click the **Complete Service** button to update the unit's status to Available.

Spectra Knee CPM Unavailable - In Service Create Patient Agreement View All Activity Edit Product

This product is **Active** and currently **In Service**

Product is being serviced
Add service notes or complete service below

There are no service notes

Add Service Note

Complete Service

Details

Asset Tag


Serial Number
12345

Product
Spectra Knee CPM

Manufacturer
DIQ - Donjoy

OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

While in the product is in service, select **Add Service Note** to record details.




The **Serviced** and **Completed Service** statuses will show on the **All Activity** page and include the User and Dates it was updated




11. Inactivating a Rental Product


If a rental unit is no longer available for rental, it can be set to Inactive. This will remove the unit from the Rental Product list page. Click the **Edit Product** button.



Shoulder CPM Available




On the **Edit** page, **uncheck the Active?** box, add an **Inactive Reason** (mandatory field) and click **Save Changes** button




OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

A **Success** message will display at the top of the Rental Product page, the unit will show it's Inactive and the reason the unit is inactive will be displayed



From the Rental Products list page, setting the filter to Inactive will display all Inactive units



12. Additional Information


Historical Rental Data



All rental transaction details (status update) that were completed in the system prior to these enhancements can be found within the Rental Product page by clicking the **View All Activity** button.

The **All Activity** page will show all historical details and each status change completed prior to the new enhancements.

New rental cycle details will continue to be added to the All Activity page, but also show as individual cycles on the main Rental Product page for better clarity




OVERVIEW OF RENTAL PROCESS IN MOTIONMD®

Adding Notes to a Rental Product

Notes specific to that unit can be added when creating a new Rental Product or by clicking the **Edit Product** button on the Rental Product page.



On **Edit** page, under the **Notes** section, type in and specific details Click **Save Changes** button



Added notes will show under the Details section on the Rental Product page and in line on the Rental Product List page.

